

UPCERA



C41

4-Axis Wet Milling Machine

Reliable Materials Expert




C41

4-Axis Wet milling Machine

- Compact and ergonomic design
- Built-in automatic circulating water system
- Rapid and accurate production with added efficiency
- Versatile indications up to 40 mm block
- Ideal for lab and clinic for same-day restorations

Materials



Device Parameters

Specification	Parameters
Dimensions	45*35*42 cm
Weight	45 KG
Voltage	220V/50Hz ; 110V/60Hz(Optional)
Power	1 KW
Compatible Compressed Air	0.42-0.8 MPa
Spindle Speed	80,000 rpm
Precision	< 0.01 mm
Number of Bur	5
Cooling Type	Water cooling
Tilt Angle	A axis: $\pm 360^{\circ}$
Materials	Glass ceramic, Composite Resin Block/Titanium block(Optional)
Application	Veneer, full contour crown, inlay, onlay, pre-milled abutment
Milling Speed	PMMA (12 minutes/pc), Glass-ceramic (10~30 minutes/pc) Pre-milled abutment (40 minutes/pc)